


CIRCULAR CUTTERS

To cut accurate samples of **CORRUGATED BOARD** for FCT and Basis Weight determination



Cut depth to be adjustable easily for different sample thickness

The blade remains totally covered by the security chassis

With three fix pins to the sample, placed in the base, to avoid displacement of the equipment during the cut operation

Sample cut up to 10 mm thickness

Four models available:

Code 5351 – **100 cm² area**

Code 5352 – **50 cm² area**

Code 6344 – **10 inch² area**

Code 6345 – **5 inch² area**

DIMENSIONS AND WEIGHT:

Dimensions: 150 x 150 x 190 mm (Ancho x Fondo x Alto)
Bos for Transport: 300 x 300 x 350 mm (Ancho x Fondo x Alto)
Weight Net/Gross: 1,5 Kg / 5 Kg

DELIVERY CONTENT:

> Circular Sample Cutter

