

DUBLE BLADE CUTTERS

Designed to paper samples, quickly and with just one operation of a certain width and 260mm length. With a fixed blade and two movable cutter blades, made of special steel.

According to standard: ISO 1924 – SCAN P38 – TAPPI...



Code: T-3159 - PRECISION DOUBLE BLADE CUTTER

- Cutting wide to be chosen between 10 - 12,7 - 15 - 20 or 25,4 mm
- Size :: 455 x 285 x 270 mm (width x depth x height)
- Net weight : 8 Kg

Code: T-3161 - PRECISION DOUBLE BLADE CUTTER

- Cutting wide : 50 mm
- Size :: 455 x 285 x 270 mm (width x depth x height)
- Net weight : 10 Kg

Before cutting, press the lever which meets the fixed cutter blade at the guide end. This clamps and tightens the sample and therefore ensures that high accuracy is obtained.

Equipment is provided with methacrylate security cover protection. This avoids the operator inserting his hands between the blades and the cutting base during cutting.