


PAPER SHEET PRESS

Press designed for the process of press-drying of the wet sheet of paper, obtained in the FH.01 sheet former.



HYDRAULIC PRESS OF MANUAL OPERATION

Superior circular plate and oscillating, to facilitate the placemen of the samples

- Capacity to 5 Kg/cm²
- Manual operation pump
- For 215 mm paper samples of diameter
- Para hojas de 215 mm de diámetro
- Gauge of control of the pressure

CONNECTIONS:

DIMENSIONS AND WEIGHT

Dimensions: 400 x 560 x 460 mm (W x D x H)

Box for Transport :550 x 650 x 550 mm (W x D x H)

Weight Net/Gross 65 Kg / 87 Kg

DELIVERY CONTENT:

> Hydraulic Manual Press

