




“Data in Real Time...”



LYNX Test Programme CMT - RCT - CCT - ECT - FCT - PAT Compression Tests.



ECT



FCT



SAMPLES CRUSH TESTER



It can work autonomously or Integrated in **LYNX Management Module** (one-to-one connection) or in the **LYNX Integral Management System Plus or Pro** (Multiple test equipment connection in LAN) with other parameters generated in the Testing Laboratory

- Quickness in carrying out tests and in obtaining results.
- No human mistakes.
- Traceability according to ISO 9000

Through this Software Test Module and using a Samples Crush Tester with its correspondent devices, **CMT – RCT – CCT - ECT – FCT – PAT** Compression tests in Liner Paper-Medium and Corrugated Board samples can be carried out. The data are displayed numerically and graphically on the PC monitor

- It marks **Maximum – Minimum – Medium Value** and **Standard Deflection**
- **Up to 100 Tests or Measurement Capacity** for each report
- **SAVE – PRINT and COPY to OFFICE CLIPBOARD FUNCTIONS**

OPTIONAL : Code T-6879 - LYNX MANAGEMENT MODULE with Basic Statistics

If you need to manage the data generated, we recommend the purchase of this module which allows further choice of language, reporting, to have a library of units, levels of password, introduction of minimum values, maximum and optimal to exploit data with statistics, charts bars, Bells of GAUSS, Comparison tolerances ..., export data to Word - Excel ..., PDF generation ...

