

PRECISION BALANCE


To make precision weighing in the laboratory



COMPATIBLE



"Data in Real Time..."



- Time measurement (typical): 2 seconds
- Adaptation to environmental conditions and Implementation
- 4 scales optimized filters
- Net Weight (approx.): 4 lbs.
- 10 units selectable weight:
G; ct; lb; oz; NG; dwt; tlc; mom; pcs
- Connecting to the network via AC adapter 220 V
- Integrated RS-232
- Manufactured in accordance with ISO 9001
- CE Mark

FEATURES:

Code	Model	Capacity (g)	Resolution (g)	Repeativity (g)	Plate size (mm)
6808	S-603	200	0,001	0,001	114 Ø