


To make precision weighing in the laboratory



- ❑ Time measurement (typical): 2 seconds
- ❑ Adaptation to environmental conditions and Implementation
- ❑ 4 scales optimized filters
- ❑ Net Weight (approx.): 4 lbs.
- ❑ 10 units selectable weight:
G; ct; lb; ozt; oz; NG; dwt; tlc; mom; pcs
- ❑ Connecting to the network via AC adapter 220 V
- ❑ Integrated RS-232
- ❑ Manufactured in accordance with ISO 9001
- ❑ CE Mark

FEATURES:

Code	Model	Capacity (g)	Resolution (g)	Repeativity (g)	Plate size (mm)
6808	S-603	200	0,001	0,001	114 Ø