

BASIS WEIGHT BALANCE

To determine the Basis Weight (grammage) in Paper, Board...



Precision electronic balance, dual display of direct reading:

*** Real Weighing:**

- Capacity: 150 g.
- Reading resolution: 0,01 g.

*** Basis Weight:**

- Capacity: 15.000 g / m²
- Reading resolution: 1 g / m²

The balance is provided with a template to cut samples of 100 x 100 mm. Through the keyboard can calculate the linear weight / m of roll. Based on this data and knowing the weight of the roll, determining the length in meters of the same. (Not available RS-232 output)